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INGLESE
What are THERMOTHERAPY products used for?
Thermotherapy is a therapeutic method which has been used for centuries as a remedy to improve blood circulation and thus, the supply of oxygen to muscles, and to reduce joint stiffness. In painful conditions resulting from minor trauma or contusions, thermotherapy can relieve pain by releasing heat onto the affected area. The medical device is also indicated for relaxing muscle tension caused by fatigue or simply cold-related cramps.
What is in the self-heating devices (white sachet)?
The self-heating devices contain a mixture of natural elements, such as iron, active carbon, salt and water. Once the hermetically sealed pack is opened, the ingredients react with oxygen in the air. The oxidation of iron generates therapeutic heat.
How do I activate the self-heating devices?
Take the device out of the special individual pack immediately before use. To 'activate' the device, hold it in your hand and wait for a few minutes. If you are using a patch, remove the protective film and apply the adhesive area on the part of the body where you feel pain or discomfort. The device will gradually heat up until it reaches the therapeutic heat (about 40°C). Do not leave the device on the table while waiting for it to heat up. Do not attempt to heat the patch in any other way.
How long does it take the product to activate?
The self-heating devices have been developed to gradually reach the therapeutic temperature to avoid thermal shock on the skin and deliver constant heat. Activate them by following the instructions outlined in the previous question or on the packaging. The self-heating devices will reach the therapeutic temperature in about 30 minutes. This time may vary depending on climate conditions, body temperature and correct activation method.

How long does the heat last?
The self-heating devices retain their heat for 8 hours.
What should I do if the patch comes off?
The natural zinc oxide-based adhesive is delicate and ideal for pain-free removal of the patch. If the patch comes off prematurely or does
not adhere well to the skin, due to a specific skin condition, you can put it back on over a thin garment to continue the treatment.
Can the self-heating devices be reused?
The self-heating devices are for single use. Once the thermal reaction has been activated, the device cannot be reused.
Can the neck and back bands be reused?
Yes, the neck and back bands can be reused. The self-heating device (refill) must be replaced for practical and efficient reuse.
When should the self-heating devices not be used for thermotherapy?
It is important to following the WARNINGS indicated on each individual pack. In particular, pay special attention to all the cases in which the self-heating patch MUST NOT be used.
Can the heat produced by the selfheating devices cause a burn?
The self-heating devices must be used in compliance with the information indicated in the WARNINGS.
Incorrect use of the patch could cause a burn. Check skin frequently during use. Remove the patch if excessive heat or itching is perceived, or if even mild signs of burning (e.g. redness) appear on the treated area.

Can I use the self-heating devices during the night?
The self-heating devices must not be used while sleeping. This is because you should avoid use in conditions that could place pressure on the sachets. You also have a reduced perception of heat and cannot monitor the correct use of the device.
Do the self-heating devices have an expiry date?
All the self-heating devices expire after 3 years. They should not be used after the expiry date.
Do the self-heating devices contain medication?
The self-heating devices are medical devices that do not contain medication.
Can I use the self-heating patch over clothing?
The patch can also be applied over an item of underwear or a thin garment. This use is recommended for the elderly as the risk of burns increases as you get older. Generally speaking, a lightweight garment acts as a protective barrier between the skin and the patch, while still allowing the diffusion of the therapeutic heat.
Does using the medical device with the band improve tolerance for people with sensitive skin and the elderly?
Yes, the medical device applied with the band is generally well tolerated, especially with particularly sensitive people or the elderly. The band allows a gentler, more comfortable application, reducing possible irritation due to direct contact of the device with the skin.
Who produces the self-heating devices?
These devices are produced in Italy by Diva International Via dell'Industria, 7 – 06038 Spello (PG).

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